




TA-04 Mouse Anti-HA mAb

Immunogen	Synthetic peptide from influenza hemaglutinin (YPYDVPDYA) coupled to KLH
Product name	Mouse anti-HA monoclonal antibody
Clone	CB051
Isotype	IgG1
Amount	0.1 ml
Reactivity	HA tag
Storage/Stability	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Avoid repeated freeze-thaws.
Application	WB
Suggested dilutions	1:1000-1:2000

Experimental data



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or HA tagged LGALS3 cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HA.

Notes: This product is to be used for laboratory only. Not for diagnostic use.



