

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA302089

FGFR4 Mouse Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommend Dilution:	WB: 1:2000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1, kappa
Clonality:	Monoclonal
Immunogen:	Purified recombinant FGFR4 fusion protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	0.25 mg/ml
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Predicted Protein Size:	87954 Da
Gene Name:	fibroblast growth factor receptor 4
Database Link:	<u>NP_002002 Entrez Gene 2264 Human</u>
Background:	FGFR4 is a member of the fibroblast growth factor receptor family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein would consist of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member preferentially binds acidic fibroblast growth factor and, although its specific function is unknown, it is overexpressed in gynecological tumor samples, suggesting a role in breast and ovarian tumorigenesis.
Synonyms:	CD334; JTK2; TKF
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Endocytosis, MAPK signaling pathway, Regulation of actin cytoskeleton



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western blot analysis of anti-FGFR4 Monoclonal Antibody (Cat.TA302089) by FGFR4 recombinant protein (Fragment 42KD). FGFR4 protein (arrow) was detected using the purified Mab.

Western blot analysis of FGFR4 (arrow) using mouse monoclonal FGFR4 antibody. 293 cell lysates (2 ug/lane) either nontransfected (Lane 1) or transiently transfected with the FGFR4 gene (Lane 2) (Origene Technologies)

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US