

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA309469

Rhodopsin (RHO) Mouse Monoclonal Antibody [Clone ID: 4D2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	4D2
Applications:	IF
Recommend Dilution:	ICC, IHC, IP, WB, ELISA
Reactivity:	Bovine, Human, Mouse, Porcine, Rat
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Bovine Rhodopsin
Formulation:	TBS, 50% glycerol, 0.09% sodium azide
Concentration:	250ug/mL
Gene Name:	rhodopsin
Database Link:	<u>NP_000530 Entrez Gene 24717 RatEntrez Gene 212541 MouseEntrez Gene 6010 Human</u>
Background:	Retinitis pigmentosa is an inherited progressive disease which is a major cause of blindness in western communities. It can be inherited as an autosomal dominant, autosomal recessive, or X-linked recessive disorder. In the autosomal dominant form, which comprises about 25% of total cases, approximately 30% of families have mutations in the gene encoding the rod photoreceptor-specific protein rhodopsin. This is the transmembrane protein which, when photoexcited, initiates the visual transduction cascade. Defects in this gene are also one of the causes of congenital stationary night blindness. [provided by RefSeq]
Synonyms:	CSNBAD1; OPN2; RP4
Protein Families:	Druggable Genome, Transmembrane



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Immunofluorescence of Rho using a dilution of Rhodopsin antibody

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US