

Product datasheet for TA309560

MRP1 (ABCC1) Mouse Monoclonal Antibody [Clone ID: IU5C1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	IU5C1
Applications:	IF, WB
Recommend Dilution:	WB, IF
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	A recombinant peptide corresponding to residues 1-33 (N-terminal: SMALRGFCSADGSDPLWDWNVTWNTSNPDFTKCF) of the human protein. [Swiss-Prot# P33527]
Formulation:	0.1% Sodium Azide
Purification:	Ascites
Gene Name:	ATP binding cassette subfamily C member 1
Database Link:	NP_004987 Entrez Gene 17250 Mouse Entrez Gene 4363 Human
Synonyms:	ABC29; ABCC; GS-X; MRP; MRP1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ABC transporters




[View online »](#)

Product images:



Western Blot: MRP1 Antibody (IU5C1) -Detection of MRP1 in 10 ug of 293 transfected (lane 2) and empty vector (lane 1) lysate.



Immunofluorescence: MRP1 Antibody (IU5C1) - Immunofluorescent staining of MRP1 using anti-MRP1 in HEK293 cells with MRP1 expression.