

Product datasheet for TA309593

BRIP1 Mouse Monoclonal Antibody [Clone ID: pp15-IB4]

Product data:

Product Type: Primary Antibodies

Clone Name: pp15-IB4

Applications: WB

Recommend Dilution: WB

Reactivity: Human Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: Partial Brip1 protein peptide.

Formulation: 0.02% Sodium Azide

Concentration: 7.4mg/ml **Purification:** Ascites

Gene Name: BRCA1 interacting protein C-terminal helicase 1

Database Link: NP 114432 Entrez Gene 83990 Human

Synonyms: BACH1; FANCJ; OF


Protein Families: Druggable Genome

Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn





Product images:



BRIP1 detected in 293 cell lysate using anti-BRIP1.