

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA309622

alpha 5 Defensin (DEFA5) Mouse Monoclonal Antibody [Clone ID: 8C8]

Product data:

Product Type: Primary Antibodies

Clone Name: 8C8
Applications: IHC

Recommend Dilution:WB, IF, IHCReactivity:HumanHost:MouseIsotype:IgG2b

Clonality: Monoclonal

Immunogen: Full-length proHD5 protein. [Swiss-Prot# Q01523] The propeptide is the full-length protein. In

the mature form, it is cleaved at residue 63 (alanine).

Formulation: Tris-glycine, 150mM NaCl and 0.05% Sodium Azide

Concentration: 2.35mg/ml

Purification: Protein G purified

Gene Name: defensin alpha 5


Database Link: NP 066290 Entrez Gene 1670 Human

Synonyms: DEF5; HD-5

Protein Families: Secreted Protein



Product images:



Immunohistochemistry: Defensin alpha 5 Antibody (8C8) - Human skeletal muscle myocytes showing moderate nuclear and faint cytoplasmic staining.