

Product datasheet for **TA309692**


HIF1AN Mouse Monoclonal Antibody [Clone ID: 162c]

Product data:


Product Type:	Primary Antibodies
Clone Name:	162c
Applications:	IHC, WB
Recommend Dilution:	WB, IHC, IF, ICC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human Factor Inhibiting HIF-1 protein expressed in E.coli BL21 (DE3) cells. [Swiss-Prot# Q9NWT6]
Formulation:	Tris-glycine, 150mM NaCl pH 7.5 and 0.05% Sodium Azide
Concentration:	1 mg/ml
Purification:	Protein G purified
Gene Name:	hypoxia inducible factor 1 alpha subunit inhibitor
Database Link:	NP_060372 Entrez Gene 55662 Human
Synonyms:	FIH1
Protein Families:	Druggable Genome



[View online »](#)

Product images:

Western Blot: Factor Inhibiting HIF-1 Antibody (FIH 162c) - Western blot analysis in A431 cell lysates.



Immunohistochemistry: Factor Inhibiting HIF-1 Antibody (162c) - IHC analysis of FIH in human renal cancer using DAB with hematoxylin counterstain.