

Product datasheet for **TA309694**

HIP1 Mouse Monoclonal Antibody [Clone ID: 1B11]


Product data:

Product Type:	Primary Antibodies
Clone Name:	1B11
Applications:	WB
Recommend Dilution:	WB, IF
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Human 3'-HIP1 (~65 kDa)
Formulation:	0.1% Sodium Azide
Purification:	Ascites
Gene Name:	huntingtin interacting protein 1
Database Link:	NP_005329 Entrez Gene 215114 Mouse Entrez Gene 3092 Human
Synonyms:	HIP-I; ILWEQ; SHON; SHONbeta; SHONgamma
Protein Families:	Druggable Genome
Protein Pathways:	Huntington's disease



[View online »](#)

Product images:



Western Blot: HIP1 Antibody (1B11) - Western Blot on HeLa whole cell extracts.