

Product datasheet for TA309695

HIP1 Mouse Monoclonal Antibody [Clone ID: 4B10]

Product data:

Product Type: Primary Antibodies

Clone Name: 4B10
Applications: WB
Recommend Dilution: WB

Reactivity: Human Host: Mouse

Isotype: IgG1, kappa **Clonality:** Monoclonal

Immunogen: Human 3'-HIP1 (~65 kDa)

Formulation: 0.1% Sodium Azide


Purification: Ascites

Gene Name: huntingtin interacting protein 1

Database Link: NP 005329 Entrez Gene 3092 Human

Synonyms: HIP-I; ILWEQ; SHON; SHONbeta; SHONgamma

Protein Families: Druggable Genome
Protein Pathways: Huntington's disease




9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn





Product images:



Western Blot: HIP1 Antibody (4B10) - Western Blot on HeLa whole cell extracts.