

Product datasheet for TA309826

RAD51C Mouse Monoclonal Antibody [Clone ID: 2H11/6]

Product data:

Product Type: Primary Antibodies

Clone Name: 2H11/6

Applications: WB

Recommend Dilution: WB

Reactivity: Human, Mouse

Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: His-tagged human Rad51C, over-expressed in E. coli.

Formulation: PBS (physiological pH) and 0.02% Sodium Azide

Purification: IgG purified

Gene Name: RAD51 paralog C

Database Link: NP 002867 Entrez Gene 114714 MouseEntrez Gene 5889 Human

Synonyms: BROVCA3; FANCO; R51H3; RAD51L2

Protein Families: Druggable Genome

Protein Pathways: Homologous recombination

OriGene Technologies, Inc.


Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com

EU: info-de@origene.com CN: techsupport@origene.cn

9620 Medical Center Drive, Ste 200



Product images:



Western Blot: Rad51C Antibody (2H11/6) - Rad51C detected in HEK293 lysate using anti-Rad51C. Photo courtesy of B.T. Bennett & K. Knight, University of Massachusetts Medical School.