


Product datasheet for TA309941

GNRHR Mouse Monoclonal Antibody [Clone ID: AT2.G7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT2.G7
Applications:	WB
Recommend Dilution:	WB: 1:1000, IHC: 1:200
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A mixture of peptides from different regions of GnRH receptor protein.
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	1 mg/ml
Gene Name:	gonadotropin releasing hormone receptor
Database Link:	NP_000397 Entrez Gene 2798 Human
Synonyms:	GNRHR1; GRHR; HH7; LHRHR; LRHR
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	GnRH signaling pathway, Neuroactive ligand-receptor interaction

Product images:



Western blot analysis of the molecular weight and levels of gonadotropin-releasing hormone (GnRH) receptor protein. Lanes 1 and 2: tissue extracts from different stages of ovarian cancers. Lane M: protein markers, molecular weight in kDa.



[View online »](#)