

Product datasheet for TA309942

ABCA1 Mouse Monoclonal Antibody [Clone ID: AB.H10]

Product data:

Product Type: Primary Antibodies

Clone Name: AB.H10

Applications: WB

Recommend Dilution: WB: 1:1000, ELISA, IP, IHC

Reactivity: Chicken, Human, Mouse

Host: Mouse Isotype: IgG

Clonality: Monoclonal

Immunogen: Peptides amino acids 1800-2260 (C-Terminal end)

Formulation: PBS, pH 7.4 with 0.05% sodium azide.

Concentration: 1 mg/ml

Gene Name: ATP binding cassette subfamily A member 1


Database Link: NP 005493 Entrez Gene 11303 MouseEntrez Gene 19 Human

Synonyms: ABC-1; ABC1; CERP; HDLDT1; TGD

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: ABC transporters

Product images:



Western blot analysis of ABCA1 induction by cholesterol (+) in human fibroblast cells from different patients (1, 2, 3, 4). Simultaneous blotting with anti-actin antibody was used for protein loading control.

OriGene Technologies, Inc.

Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

9620 Medical Center Drive, Ste 200

