

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA309943

CD133 (PROM1) Mouse Monoclonal Antibody [Clone ID: EMK08]

Product data:

Product Type: Primary Antibodies

Clone Name: EMK08

Applications: WB

Recommend Dilution: WB 1:500
Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Homo sapiens prominin 1 (PROM1), transcript variant 1 whole protein (Accession ref:

NP 006008)

Formulation: PBS, pH 7.4 with less than or equal to 0.09% sodium azide.

Purification: Affinity purified

Gene Name: prominin 1


Database Link: NP 006008 Entrez Gene 8842 Human

Synonyms: AC133; CD133; CORD12; MCDR2; MSTP061; PROML1; RP41; STGD4

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane



Product images:



Western blot analysis of HCT116 WT and HCT116 p53-/- cell lines by anti-CD133 antibody (1:500 dilution on 100ug of cell lysates). Left: HCT116 - Human colorectal epithelial carcinoma (expressing CD133); Right: p53 knock out line - no expression of CD133.