

Product datasheet for TA317480


UBA6 Mouse Monoclonal Antibody [Clone ID: 1D11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1D11
Applications:	IHC
Recommend Dilution:	IHC-P (5 µg/ml)
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	UBE1L2 / UBE1L2 antibody was raised against fLJ10808 (NP_060697, 962 a.a. ~ 1052 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Formulation:	PBS, pH 7.2
Concentration:	0.5mg/ml
Purification:	Protein A Column
Gene Name:	ubiquitin like modifier activating enzyme 6
Database Link:	NP_060697 Entrez Gene 55236 Human
Synonyms:	E1-L2; MOP-4; UBE1L2
Note:	Specific for Human UBA6
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)

Product images:

Anti-UBA6 antibody IHC of human testis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.