

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA319910

DC SIGN (CD209) Mouse Monoclonal Antibody [Clone ID: 8B6]

Product data:

Product Type: Primary Antibodies

Clone Name: 8B6

Applications: IHC, WB

Recommend Dilution: WB: 1 - 2 mg/mL, ICC: 5 - 10 mg/mL

Reactivity: Human Host: Mouse

Clonality: Monoclonal

Immunogen: A recombinant His-tagged protein fragment corresponding to the extracellular region of

human DC-SIGN.

Formulation: DC-SIGN Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.

Concentration: 1ug/ul

Purification: DC-SIGN Monoclonal Antibody is immunoaffinity chromotography purified IgG.

Gene Name: CD209 molecule

Database Link: NP 001138365 Entrez Gene 30835 Human

Background: DC-SIGN Monoclonal Antibody: Dendritic cells (DCs) that control immune responses were

recently found to capture and transport HIV from the mucosal area to remote lymph nodes, where DCs hand over HIV to CD4+ T lymphocytes. DCs also amplify the amount of virus and extend the duration of viral infectivity. Multiple strains of HIV-1, HIV-2 and SIV bind to DCs via DC-SIGN. ICAM-3 is the natural ligand for DC-SIGN. A DC-SIGN homologue (termed DC-SIGNR, L-SIGN, and DC-SIGN2) was identified recently. DC-SIGN forms a novel gene family with DC-SIGNR and many alternatively spliced isoforms of DC-SIGN and DC-SIGNR are known to exisit. The expression of DC-SIGN was found in mucosal tissues including placenta, small intestine,

and rectum.


Synonyms: CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1

Protein Families: Druggable Genome






Product images:



Western blot detection of DC-SIGN fusion protein in human uterus tissue lysate at (A) 1 and (B) 2 ug/mL.



Immunohistochemistry of DC-SIGN in lymph node tissue with DC-SIGN antibody at 10 ug/mL.