

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA325138

Creatine kinase B type (CKB) Mouse Monoclonal Antibody [Clone ID: 893CT29.1.1]

Product data:

Product Type: Primary Antibodies

Clone Name: 893CT29.1.1

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Formulation: Purified His-tagged CKB protein was used to produced this monoclonal antibody. Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: 0.37 mg/ml

Purification: This antibody is purified through a protein G column, followed by dialysis against PBS.

Predicted Protein Size: 42644 Da

Gene Name: creatine kinase B

Database Link: NP 001814 Entrez Gene 12709 MouseEntrez Gene 1152 Human

Synonyms: B-CK; BCK; CKBB; HEL-211; HEL-S-29


Protein Families: Druggable Genome

Protein Pathways: Arginine and proline metabolism, Metabolic pathways





Product images:



Creatine Kinase BB (CKB) Antibody (Center) (Cat. #TA325138) western blot analysis in MDA-MB453, 293, Y79 cell line, mouse stomach and brain tissue lysates (35ug/lane). This demonstrates the CKB antibody detected the CKB protein (arrow).