

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA326443

Ankyrin brain (ANK2) Mouse Monoclonal Antibody [Clone ID: S105-13]


Product data:

Product Type:	Primary Antibodies
Clone Name:	S105-13
Applications:	WB
Recommend Dilution:	IHC: 1:50-1:200, ICC: 1:50-1:200
Reactivity:	Human, Mouse, Rat
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Synthetic peptide amino acids 203-496 of human Ankyrin-B
Formulation:	PBS pH7.4, 50% glycerol, 0.09% sodium azide
Concentration:	1 mg/ml
Purification:	Protein G Purified
Gene Name:	ankyrin 2, neuronal
Database Link:	<u>NP_001120965 Entrez Gene 109676 MouseEntrez Gene 362036 RatEntrez Gene 287 Human</u>
Background:	Ankyrins are a family of adaptor proteins that mediate the attachment of integral membrane proteins to thespectrin-actin based membrane skeleton . Ankyrins have binding sites for the beta subunit of spectrin and at least 12 families of integral membrane proteins. This linkage is required to maintain the integrity of theplasma membranes and to anchor specific ion channels, ion exchangers and ion transporters in the plasma membranes.
Synonyms:	ANK-2; brank-2; LQT4
Note:	Detects ~200kDa. No cross reactivity against Ankyrin-G
Protein Families:	Druggable Genome



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western blot analysis of Ankyrin-B in rat brain membrane tissues, using a 1:1000 dilution of the antibody

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US