

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA326482

# Kv1.2 (KCNA2) Mouse Monoclonal Antibody [Clone ID: S14-16]


## **Product data:**

Product Type:	Primary Antibodies
Clone Name:	S14-16
Applications:	IHC
<b>Recommend Dilution:</b>	WB: 1ug/ml, IHC: 0.1-1ug/ml, IF: 1-10ug/ml
Reactivity:	Human, Mouse, Rat, Xenopus, Zebrafish
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Fusion protein amino acids 428-499 (cytoplasmic C-terminus) of rat Kv1.2
Formulation:	PBS pH7.4, 50% glycerol, 0.09% sodium azide
Concentration:	1 mg/ml
Purification:	Protein G Purified
Gene Name:	potassium voltage-gated channel subfamily A member 2
Database Link:	<u>NP_004965 Entrez Gene 16490 MouseEntrez Gene 25468 RatEntrez Gene 3737 Human</u>
Background:	Voltage gated channels are tetrameters composed of four alpha-subunits arranged around a central pore. Each alphasubunit consists of six transmembrane segments with cytoplasmic NH2 and COOH-termini. Members of the KV1- KV4 subfamilies generate functional K+ channels in a homotetrameric configuration.
Synonyms:	EIEE32; HBK5; HK4; HUKIV; KV1.2; MK2; NGK1; RBK2
Note:	Detects ~80kDa. No cross-reactivity against Kv1.1, 1.3, 1.4 or 1.6 expressed in transfected cells.
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

## **Product images:**



IHC analysis of Kv1.2 in frozen sections of mouse brain extract using the antibody

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US