

Product datasheet for TA326493

MDC1 Mouse Monoclonal Antibody [Clone ID: P2B11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	P2B11
Recommend Dilution:	WB: 1:2000
Reactivity:	Bovine, Chimpanzee, Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	GST-tagged recombinant protein corresponding to mouse MDC1 ar and around the N-terminus
Formulation:	PBS pH7.4, 50% glycerol, 0.09% sodium azide
Concentration:	1 mg/ml
Purification:	Protein G Purified
Gene Name:	mediator of DNA damage checkpoint 1
Database Link:	NP_055456 Entrez Gene 240087 MouseEntrez Gene 9656 Human
Background:	MDC1, mediator of DNA damage checkpoint protein 1, plays a role in checkpoint mediated cell cycle arrest in response to DNA damage, within S phase and G2/M. It is also thought to act as a scaffold protein during recruitment of DNA repair and signal transduction proteins to discrete foci of DNA damage that are marked by phosphorylation of histone H2A.X on S139.
Synonyms:	NFBD1
Note:	Detects ~184kDa. This antibody recognizes MDC1 at and around the N-terminus.



[View online »](#)