

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA327648

CD61 (ITGB3) Mouse Monoclonal Antibody [Clone ID: 2F2]

Product data:

Product Type: Primary Antibodies

Clone Name: 2F2
Applications: IHC

Recommend Dilution: IHC: 1:25 - 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Formulation: The antibody is diluted in tris buffered saline, pH 7.3-7.7, with 1% BSA and <1% sodium azide.

Purification: Affinity purification

Gene Name: integrin subunit beta 3

Database Link: NP 000203 Entrez Gene 3690 Human

Synonyms: BDPLT2; BDPLT16; CD61; GP3A; GPIlla; GT

Note: CD61 is the human integrin beta chain beta 3 protein (ITGB3). Integrins are integral cell-

surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-

surface mediated signaling. The integrin beta 3 chain of the vitronectin receptor and GPIIb/IIIa complex is a 90-110 kDa glycoprotein polypeptide which is expressed on platelets, megakaryocytes, macrophages, osteoclasts, and synovial lining cells. Integrin alpha-IIb/beta-3 recognizes the sequence H-H-L-G-G-G-A-K-Q-A-G-D-V in fibrinogen gamma chain. Following activation, integrin alpha-IIb/beta-3 brings about platelet/platelet interaction through binding of soluble fibrinogen. This step leads to rapid platelet aggregation which physically plugs ruptured endothelial surface. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions. This antibody is

of obviously abnormal forms and micromegakaryocytes in myelodysplastic syndrome, acute myeloid leukemia with multilineage dysplasia, acute megakaryoblastic leukemia, and

useful in evaluating the number of megakaryocytes, size, nuclear lobation, and the presence

myeloproliferative neoplasms.





CD61 (ITGB3) Mouse Monoclonal Antibody [Clone ID: 2F2] - TA327648


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways: Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-

receptor interaction, Focal adhesion, Hematopoietic cell lineage, Hypertrophic

cardiomyopathy (HCM), Regulation of actin cytoskeleton

Product images:



Immunohistochemistry staining of Paraffin Bone marrow tissue by CD61 antibody (dilution: 1:25 - 1:100; visualization of staining: Cytoplasmic)