

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA327925

LAT Mouse Monoclonal Antibody [Clone ID: 11B.12]

Product data:

Product Type: Primary Antibodies

Clone Name: 11B.12 Applications: WB

Recommend Dilution: WB, IP

Reactivity: Human, Mouse

Host: Mouse

Isotype: IgG1, kappa **Clonality:** Monoclonal

Immunogen: GST fusion protein corresponding to residues 31-233 of human LAT, cytoplasmic domain

Formulation: This antibody is provided in phosphate-buffered solution, pH 7.2, containing 0.09% sodium

azide. Final antibody concentration is 0.5 mg/ml.

Concentration: 0.5 mg/ml

Purification: The antibody was purified by affinity chromatography.

Predicted Protein Size: 28 kD but expected at 36-38 kD Gene Name: linker for activation of T-cells

Database Link: NP 001014987 Entrez Gene 16797 MouseEntrez Gene 27040 Human

Background: LAT is a transmembrane molecule that mediates signaling following the engagement of a

number of cell surface receptors including TCR during mature T cell activation, pre-TCR during T cell development, Fc RIII on NK cells, and FcɛRI on mast cells.LAT is multiply phosphorylated on tyrosines by ZAP-70 and Syk, which leads to downstream effects that include release of

intracellular calcium, PKC activation, MAPK activation and cytoskeletal

reorganization.Palmitoylation of LAT is required for its localization to lipid rafts and efficient

phosphorylation.

Synonyms: LAT1; pp36

Protein Families: Druggable Genome, Transmembrane


Protein Pathways: Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Natural killer cell

mediated cytotoxicity, T cell receptor signaling pathway





Product images:



Western blot analysis of extracts from Jurkat cells (lane 1) and mouse thymocytes (lane 2) using anti-LAT, clone 11B.12.