

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA328026

# ALIX (PDCD6IP) Mouse Monoclonal Antibody [Clone ID: 3A9]

## **Product data:**


Product Type:	Primary Antibodies
Clone Name:	3A9
Applications:	WB
<b>Recommend Dilution:</b>	WB, IP, IF
Reactivity:	Frog, Human, Mouse
Host:	Mouse
lsotype:	lgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant full-length human Alix protein
Formulation:	This antibody is provided in phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Concentration:	0.5 mg/ml
Purification:	The antibody was purified by affinity chromatography.
Predicted Protein Size:	96 kD
Predicted Protein Size: Gene Name:	96 kD programmed cell death 6 interacting protein
Gene Name:	programmed cell death 6 interacting protein
Gene Name: Database Link:	programmed cell death 6 interacting protein <u>NP 037506 Entrez Gene 18571 MouseEntrez Gene 10015 Human</u> Alix (ALG-2-interacting protein X) is an adaptor protein that was first described for its capacity to bind to the calcium-binding protein ALG-2, the expression of which seemed necessary for cell death. Predicted molecular weight approximately 96 kD. Alix contains an N-terminal Bro1 domain and a C-terminal proline-rich domain (PRD). The PRD contains multiple polyproline motifs, which are potential docking sites for proteins containing SH3 domains, interactwith multiple cellular Alix-binding partners that are involvedin apoptotic induction, endosomal sorting, and endocytosis. Bro1 domain interacts with CHMP4b which is also involved in
Gene Name: Database Link: Background:	programmed cell death 6 interacting protein <u>NP 037506 Entrez Gene 18571 MouseEntrez Gene 10015 Human</u> Alix (ALG-2-interacting protein X) is an adaptor protein that was first described for its capacity to bind to the calcium-binding protein ALG-2, the expression of which seemed necessary for cell death. Predicted molecular weight approximately 96 kD. Alix contains an N-terminal Bro1 domain and a C-terminal proline-rich domain (PRD). The PRD contains multiple polyproline motifs, which are potential docking sites for proteins containing SH3 domains, interactwith multiple cellular Alix-binding partners that are involvedin apoptotic induction, endosomal sorting, and endocytosis. Bro1 domain interacts with CHMP4b which is also involved in endosomal sorting.



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

### **Product images:**



NIH3T3 whole cell extracts were resolved by electrophoresis, transferred to nitrocellulose, and probed with anti-Alix antibody (clone 3A9). Proteins were visualized using a goat anti-mouse-IgG secondary conjugated to HRP and chemiluninescence detection.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US