

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336397

IL17E (IL25) Mouse Monoclonal Antibody [Clone ID: 68C1039.2]

Product data:

Product Type: Primary Antibodies

Clone Name: 68C1039.2

Applications: IHC, WB

Recommend Dilution: WB: 2-4 ug/ml, IHC: 1:10-1:500, IHC-P: 1:10-1:500

Reactivity: Human, Mouse

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Amino acids 115-132 of human IL-17E (isoform CRA_a) were used as immunogen for this

antibody.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Predicted Protein Size: 20 kDa

Gene Name: interleukin 25

Database Link: NP 758525 Entrez Gene 140806 MouseEntrez Gene 64806 Human

Background: A number of cytokines belonging to the interleukin (IL)-17 family have been identified, termed

IL-17A-F. IL-17 is a potent proinflammatory cytokine that plays roles in a number of diseases including rheumatoid arthritis, multiple sclerosis, and promotion of tumor growth. IL-17E, C, and B are able to induce proinflammatory responses. However, they do not bind to the IL-17 receptor suggesting that an additional IL-17R related receptor may exist. Receptors for IL-17B and IL-17E have been independently isolated by Shi, et al and Lee, et al. and have been designated as EV127 (in mouse) and IL-17Rh1 (in human), respectively. IL-17E, also called IL-

25, induces activation of NF-kB pathway and like IL-17 also induces production of IL-8.

Synonyms: IL17E

Note: IHC-Paraffin reported in literature (see Terrier B, et al).


Protein Families: Druggable Genome, Secreted Protein, Transmembrane






Protein Pathways: Cytokine-cytokine receptor interaction

Product images:



Western Blot: IL17E Antibody (68C1039.2) TA336397 - Analysis of IL-17E using IL-17E monoclonal antibody. Mouse testis lysate probed with IL-17E antibody at 2 ug/ml.



Immunohistochemistry: IL17E Antibody (68C1039.2) TA336397 - Adrenal cortex, Human: Formalin-Fixed Paraffin-Embedded (FFPE)