

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA336944

Annexin A1 (ANXA1) Mouse Monoclonal Antibody [Clone ID: 2F1]

Product data:

Product Type: Primary Antibodies

Clone Name: 2F1

Applications: IHC, WB

Recommend Dilution: WB: 1:500 - 1:2000, ELISA: 1:10000, IHC: 1:200 - 1:1000, IHC-P: 1:200 - 1:1000

Reactivity: Human, Mouse, Primate

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Partial recombinant human Annexin A1 (between residues 100-300) expressed in E. coli

[UniProt P04083]

Formulation: PBS, 0.03% Sodium Azide. Store at 4C short term. Aliquot and store at -20C long term. Avoid

freeze-thaw cycles.

Concentration: 1 mg/ml

Purification: Ammonium sulfate precipitation

Gene Name: annexin A1

Database Link: NP 000691 Entrez Gene 16952 MouseEntrez Gene 301 Human

Background: Annexin I belongs to a family of Ca(2+)-dependent phospholipid binding proteins which have

a molecular weight of approximately 35,000 to 40,000 and are preferentially located on the cytosolic face of the plasma membrane. Annexin I protein has an apparent relative molecular mass of 40 kDa, with phospholipase A2 inhibitory activity. Since phospholipase A2 is required

for the biosynthesis of the potent mediators of inflammation, prostaglandins and

leukotrienes, annexin I may have potential anti-inflammatory activity. [provided by RefSeq]


Synonyms: ANX1; LPC1

Protein Families: Druggable Genome






Product images:




Western Blot: Annexin A1 Antibody (2F1) [TA336944] - WB analysis of ANXA1 in (1) Hela, (2) A549, (3) K562, (4) NIH3T3, (5) C6, and (6) COS7 cell lysate.



Western Blot: Annexin A1 Antibody (2F1) [TA336944] - WB detection of ANXA1 in human ANXA1 recombinant protein.





Immunohistochemistry: Annexin A1 Antibody (2F1) [TA336944] - IHC analysis of paraffin embedded human cervical cancer tissues with DAB staining.