

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA337067

TLR5 Mouse Monoclonal Antibody [Clone ID: 19D759.2]

Product data:

Product Type: Primary Antibodies

Clone Name: 19D759.2

Applications: FC, IHC, WB

Recommend Dilution: WB: 1-3 ug/ml, FC: 0.5-1 ug/10^6 cells, IHC-F: 10 ug/ml, IHC-P: 5 ug/ml

Reactivity: Canine, Human, Mouse

Host: Mouse

Isotype: IgG2a, kappa
Clonality: Monoclonal

Immunogen: This antibody was developed against KLH-conjugated synthetic peptide corresponding to a

portion of human TLR5 between amino acids 700-800. The antibody also reacts with mouse

TLR5.

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at -20C. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Gene Name: toll like receptor 5

Database Link: NP 003259 Entrez Gene 53791 MouseEntrez Gene 7100 Human

Synonyms: MELIOS; SLE1; SLEB1; TIL3

Note: Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM sodium

citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 min.

Protein Families: Druggable Genome, Transmembrane


Protein Pathways: Pathogenic Escherichia coli infection, Toll-like receptor signaling pathway






Product images:

Validation of NBP2-24787 anti-TLR5 mouse monoclonal antibody.




Western Blot: TLR5 Antibody (19D759.2) TA337067 - BEAS-2B cells whole cell lysate. Image from verified customer review.




Immunohistochemistry-Paraffin: TLR5 Antibody (19D759.2) TA337067 - analysis of TLR5 in formalin-fixed, paraffin-embedded mouse brain tissue using TLR5 antibody at 5 ug/ml.






Immunohistochemistry-Paraffin: TLR5 Antibody (19D759.2) TA337067 - Formalin-fixed, paraffinembedded human skin stained with TLR5 antibody (5 ug/ml), peroxidase-conjugate and DAB chromogen.



Immunohistochemistry-Paraffin: TLR5 Antibody (19D759.2) TA337067 - IHC analysis of formalinfixed paraffin-embedded tissue section of normal human skin using 5 ug/ml concentration of TLR5 antibody (clone 19D759.2). Intense membrane staining of TLR5 was



Flow Cytometry: TLR5 Antibody (19D759.2) TA337067 - Intracellular flow analysis of TLR5 in human monocytes using 0.5 ug of TLR5 antibody. Shaded histogram represents cells without antibody; green represents isotype control; red represents TLR5 antibody.