

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA337075**

## TLR3 Mouse Monoclonal Antibody [Clone ID: 40C1285.6]


### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	40C1285.6
Applications:	FC, IHC, WB
Recommend Dilution:	WB: 1-3 ug/ml, FC: 2-4 ug/ 1x10^6 cells, IF: 1:10-1:500, IHC: 1:10-1:500, IHC-P: 1:10-1:500, IP: 2 ug/ 1x10^6 cells
Reactivity:	Canine, Human, Mouse
Host:	Mouse
lsotype:	lgG1, kappa
Clonality:	Monoclonal
Immunogen:	A synthetic peptide corresponding to amino acids 55-85 of human TLR3 was used as immunogen.
Formulation:	PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at - 20C long term. Avoid freeze-thaw cycles.
Concentration:	0.5 mg/ml
Purification:	Protein G purified
Gene Name:	toll like receptor 3
Database Link:	NP 003256 Entrez Gene 142980 MouseEntrez Gene 7098 Human
Synonyms:	CD283; IIAE2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Toll-like receptor signaling pathway




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### **Product images:**




Western Blot: TLR3 Antibody (40C1285.6) TA337075 - analysis of TLR3 in lysates from A) untransfected 293, B) 293 cells transfected with human TLR3 cDNA, C) human intestine, D) placenta, E) heart and F) ovary using TLR3 antibody at 3 ug/ml. Goat anti-mou



Immunohistochemistry-Paraffin: TLR3 Antibody (40C1285.6) TA337075 - Analysis of TLR3 in (A) human gut lumen (longitudinal section, transverse region) using this antibody at 10 ug/ml (Data courtesy Dr. Elizabeth Furrie, University of Dundee) and mouse sp

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Flow Cytometry: TLR3 Antibody (40C1285.6) TA337075 - Intracellular flow analysis of TLR3 in human monocytes using 0.5 ug of the antibody. Shaded histogram represents monocytes without antibody; green represents isotype control (Novus, 20109); red repres

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US