

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA337171

HMG1 (HMGB1) Mouse Monoclonal Antibody [Clone ID: 19N12A1]

Product data:

Product Type: Primary Antibodies

Clone Name: 19N12A1

Applications: WB

Recommend Dilution: WB: 1-3ug/ml, Simple Western 1:100

Reactivity: Human, Mouse

Host: Mouse

Clonality: Monoclonal

Immunogen: Full-length recombinant human HMGB1 was used as the immunogen. full length recombinant

(inhouse)

Formulation: PBS containing 0.05% BSA, 0.05% Sodium Azide. Store at 4C short term. Aliquot and store at -

20C long term. Avoid freeze-thaw cycles.

Concentration: 0.5 mg/ml

Purification: Protein G purified

Gene Name: high mobility group box 1

Database Link: NP 002119 Entrez Gene 15289 MouseEntrez Gene 3146 Human

Background: HMGB1 (high mobility group protein B1) belongs to the HMGB family and contains two HMG

box DNA-binding domains. It is a highly conserved, ubiquitous protein present in the nuclei

and cytoplasm of nearly all cell types, and is a necessary and sufficient mediator of inflammation during sterile and infection-associated responses. HMGB1 also act as DNA nuclear binding protein that has recently been shown to be an early trigger of sterile inflammation in animal models of trauma-hemorrhage via the activation of the Toll-like receptor 4 (TLR4) and the receptor for the advanced glycation endproducts (RAGE). It is reported that the level of HMGB1 is elevated during sterile tissue injury, infection, lethal endotoxemia or sepsis, collagen-induced arthritis, and ischemia-reperfusion induced tissue

injury.

Synonyms: HMG1; HMG3; SBP-1


Protein Families: Druggable Genome, Stem cell - Pluripotency, Transcription Factors

Protein Pathways: Base excision repair






Product images:



Simple Western: HMGB1 Antibody (19N12A1) TA337171 - Simple Western lane view shows a specific band detected for HMGB1 in HeLa lysate. This experiment was performed under reducing conditions using the Wes or Sally Sue separation system 12-230kDa (or 66-



Western Blot: HMGB1 Antibody (19N12A1) TA337171 - Analysis using HMGB1 antibody. A) Full-length human HMGB1 protein, B) human HeLa cell lysate and C) mouse NIH 3T3 cell lysate were probed with HMGB1 antibody at 2 ug/ml. goat anti-mouse Ig HRP secondary