

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA346951

FYN Mouse Monoclonal Antibody [Clone ID: 4B8-E7-A8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	4B8-E7-A8
Applications:	IF, WB
Recommend Dilution:	WB: 1:500, IF: 1:50
Reactivity:	Human, Monkey, Mouse, Rat
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	The immunogen for FYN antibody: purified recombinant human Fyn protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification:	Affinity purified
Predicted Protein Size:	59 kDa
Gene Name:	FYN proto-oncogene, Src family tyrosine kinase
Database Link:	<u>NP_002028 Entrez Gene 14360 MouseEntrez Gene 25150 RatEntrez Gene 2534 Human</u>
Background:	This gene is a member of the protein-tyrosine kinase oncogene family.It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth.The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein.Alternatively spliced transcript variants encoding distinct isoforms exist. [provided by RefSeq,Jul 2008]
Synonyms:	p59-FYN; SLK; SYN
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Adherens junction, Axon guidance, Fc epsilon RI signaling pathway, Focal adhesion, Natural killer cell mediated cytotoxicity, Pathogenic Escherichia coli infection, Prion diseases, T cell receptor signaling pathway, Viral myocarditis




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



Western blot detection of Fyn in Hela, 3T3, C6, COS7, CEM, Ramos and Jurkat cell lysates using Fyn mouse mAb (1:500 diluted).Predicted band size:59KDa.Observed band size:59KDa.



Immunocytochemistry staining of HeLa cells fixed with -20? Methanol and using Fyn mouse mAb (dilution 1:50).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US