

Product datasheet for TA346964

Hexokinase 1 (HK1) Mouse Monoclonal Antibody [Clone ID: 6B8-G11-E7]

Product data:


Product Type:	Primary Antibodies
Clone Name:	6B8-G11-E7
Applications:	WB
Recommend Dilution:	WB: 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for HK1 antibody: purified recombinant human Hexokinase 1 protein fragments expressed in E.coli
Formulation:	Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification:	Affinity purified
Predicted Protein Size:	102 kDa
Gene Name:	hexokinase 1
Database Link:	NP_000179 Entrez Gene 15275 Mouse Entrez Gene 25058 Rat Entrez Gene 3098 Human
Background:	Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in this gene have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of this gene results in five transcript variants which encode different isoforms, some of which are tissue-specific. Each isoform has a distinct N-terminus; the remainder of the protein is identical among all the isoforms. A sixth transcript variant has been described, but due to the presence of several stop codons, it is not thought to encode a protein. (provided by RefSeq, Apr 2009).
Synonyms:	hexokinase; HK; HK1-ta; HK1-tb; HK1-tc; HKD; HKI; HMSNR; HXK1
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway, Metabolic pathways, Starch and sucrose metabolism, Type II diabetes mellitus

Product images:



Western blot detection of Hexokinase 1 in MCF7, Jurkat, CHO-K1, C2C12, C6, 3T3 and HeLa cell lysates using Hexokinase 1 mouse mAb (1:1000 diluted). Predicted band size:102KDa. Observed band size:102KDa.