

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA346976

IMP3 (IGF2BP3) Mouse Monoclonal Antibody [Clone ID: 1F1-E10-D11]


Product data:

Product Type:	Primary Antibodies
Clone Name:	1F1-E10-D11
Applications:	WB
Recommend Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	The immunogen for IGF2BP3 antibody: purified recombinant human IMP3 protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Purification:	Affinity purified
Predicted Protein Size:	70 kDa
Gene Name:	insulin like growth factor 2 mRNA binding protein 3
Database Link:	<u>NP_006538 Entrez Gene 10643 Human</u>
Background:	The protein encoded by this gene is primarily found in the nucleolus,where it can bind to the 5' UTR of the insulin-like growth factor II leader 3 mRNA and may repress translation of insulin-like growth factor II during late development.The encoded protein contains several KH domains,which are important in RNA binding and are known to be involved in RNA synthesis and metabolism.A pseudogene exists on chromosome 7,and there are putative pseudogenes on other chromosomes. [provided by RefSeq,Jul 2008]
Synonyms:	CT98; IMP-3; IMP3; KOC; KOC1; VICKZ3



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western blot detection of IMP3 in MCF7, SKVO3, NTERA2, Jurkat, Hela and K562 cell lysates using IMP3 mouse mAb (1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US