

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA346993

LCK Mouse Monoclonal Antibody [Clone ID: 3E5-C8-B9]

Product data:

Product Type: Primary Antibodies

Clone Name: 3E5-C8-B9

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: The immunogen for LCK antibody: purified recombinant human LCK protein fragments

expressed in E.coli.

Formulation: Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with

0.02% sodium azide, 50%, glycerol

Purification: Affinity purified

Predicted Protein Size: 58 kDa

Gene Name: LCK proto-oncogene, Src family tyrosine kinase

Database Link: NP 001036236 Entrez Gene 3932 Human

Background: This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded

protein is a key signaling molecule in the selection and maturation of developing T-cells. It contains N-terminal sites for myristylation and palmitylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-

containing and proline-rich motifs, respectively. The protein localizes to the plasma

membrane and pericentrosomal vesicles, and binds to cell surface receptors, including CD4 and CD8, and other signaling molecules. Multiple alternatively spliced variants, encoding the

same protein, have been described.

Synonyms: IMD22; LSK; p56lck; pp58lck; YT16

Protein Families: Druggable Genome, Protein Kinase, Stem cell - Pluripotency


Protein Pathways: Natural killer cell mediated cytotoxicity, Primary immunodeficiency, T cell receptor signaling

pathway





Product images:



Western blot detection of LCK in Jurkat and Ramos cell lysates and using LCK mouse mAb (1:1000 diluted). Predicted band size: 58KDa. Observed band size: 58KDa.