

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347014

OLIG2 Mouse Monoclonal Antibody [Clone ID: 7A1-B12-G4]

Product data:

Product Type: Primary Antibodies

Clone Name: 7A1-B12-G4

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: The immunogen for OLIG2 antibody: purified recombinant human OLIG2 beta protein

fragments expressed in E.coli.

Formulation: Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and

50% glycerol.

Purification: Affinity purified

Predicted Protein Size: 32 kDa

Gene Name: oligodendrocyte lineage transcription factor 2

Database Link: NP 005797 Entrez Gene 10215 Human

Background: This gene encodes a basic helix-loop-helix transcription factor which is expressed in

oligodendroglial tumors of the brain. The protein is an essential regulator of ventral neuroectodermal progenitor cell fate. The gene is involved in a chromosomal translocation t(14;21)(q11.2;q22) associated with T-cell acute lymphoblastic leukemia. Its chromosomal location is within a region of chromosome 21 which has been suggested to play a role in

learning deficits associated with Down syndrome.


Synonyms: BHLHB1; bHLHe19; OLIGO2; PRKCBP2; RACK17

Protein Families: Druggable Genome, Transcription Factors





Product images:



Western blot detection of OLIG2 in CHO-K1 cell lysate (B)and CHO-K1 transfected by EGFP-OLIG2 fragment (A)cell lysate using OLIG2 mouse mAb (1:1000 diluted).