

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347017

PI 3 Kinase p85 beta (PIK3R2) Mouse Monoclonal Antibody [Clone ID: 8D9-D5-F8]

Product data:

Product Type: Primary Antibodies

Clone Name: 8D9-D5-F8

Applications: WB

Recommend Dilution: WB: 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: The immunogen for PIK3R2 antibody: purified recombinant human PI3 Kinase p85 beta

protein fragments expressed in E.coli.

Formulation: PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Predicted Protein Size: 85 kDa

Gene Name: phosphoinositide-3-kinase regulatory subunit 2

Database Link: NP 005018 Entrez Gene 5296 Human

Background: Phosphatidylinositol 3-kinase (PI3K) is a lipid kinase that phosphorylates phosphatidylinositol

and similar compounds, creating second messengers important in growth signaling pathways. PI3K functions as a heterodimer of a regulatory and a catalytic subunit. The protein encoded by this gene is a regulatory component of PI3K. Two transcript variants, one

protein coding and the other non-protein coding, have been found for this gene.

Synonyms: MPPH; MPPH1; p85; p85-BETA; P85B

Protein Families: Druggable Genome






Protein Pathways:

Acute myeloid leukemia, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, Leukocyte transendothelial migration, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Phosphatidylinositol signaling system, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, VEGF signaling pathway

Product images:



Western blot detection of PI3 Kinase p85 beta in C6, K562, MCF7, U87MG, HCT116, 3T3 and Hela cell lysates using PI3 Kinase p85 beta mouse mAb (1:1000 diluted). Predicted band size:85KDa. Observed band size:85KDa.