

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA347044

Syntenin (SDCBP) Mouse Monoclonal Antibody [Clone ID: 3D9-G9-H4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	3D9-G9-H4
Applications:	WB
Recommend Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	The immunogen for SDCBP antibody: purified recombinant human Syntenin protein fragments expressed in E.coli.
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.02% sodium azide, 50%,glycerol
Purification:	Affinity purified
Predicted Protein Size:	32 kDa
Gene Name:	syndecan binding protein
Database Link:	<u>NP_001007068 Entrez Gene 6386 Human</u>
Background:	The protein encoded by this gene was initially identified as a molecule linking syndecan- mediated signaling to the cytoskeleton. The syntenin protein contains tandemly repeated PDZ domains that bind the cytoplasmic, C-terminal domains of a variety of transmembrane proteins. This protein may also affect cytoskeletal-membrane organization, cell adhesion, protein trafficking, and the activation of transcription factors. The protein is primarily localized to membrane-associated adherens junctions and focal adhesions but is also found at the endoplasmic reticulum and nucleus. Alternative splicing results in multiple transcript variants encoding different isoforms.
Synonyms:	MDA-9; MDA9; ST1; SYCL; TACIP18
Protein Families:	Druggable Genome, Transmembrane



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western blot detection of Syntenin in HelaNE and Hela cell lysates using Syntenin mouse mAb (1:1000 diluted).Predicted band size:32KDa.Observed band size:32KDa.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US