

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA500073**

## FOXA2 Mouse Monoclonal Antibody [Clone ID: OTI3C10]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3C10

**Applications:** FC, IF, IHC, WB

**Recommend Dilution:** WB: 1:1000; IHC: 1:50; IF: 1:50-1:100; Flow 1:1,000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 300-458 of human

FOXA2

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 48.1 kDa

**Gene Name:** forkhead box A2

Database Link: NP 710141 Entrez Gene 3170 Human

**Background:** FOXA2 is a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear

factors are transcriptional activators for liver-specific genes such as albumin and

transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. This gene has been linked to sporadic cases of maturity-onset diabetes of the young. Transcript

variants encoding different isoforms have been identified for this gene

Synonyms: HNF3B; TCF3B

Protein Families: Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription


**Factors** 

**Protein Pathways:** Maturity onset diabetes of the young







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FOXA2 ([RC211408], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FOXA2. Positive lysates [LY407005] (100ug) and [LC407005] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 1 cell line by using anti-FOXA2 monoclonal antibody at 1:200 dilution.




Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-FOXA2 monoclonal antibody (1:200).






Immunohistochemical staining of paraffinembedded Carcinoma of liver using anti-FOXA2 (TA500073) mouse monoclonal antibody (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Immunohistochemical staining of paraffinembedded lung within the normal limits using anti-FOXA2 (TA500073) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Anti-FOXA2 mouse monoclonal antibody (TA500073) immunofluorescent staining of HeLa cells transiently transfected by pCMV6-ENTRY FOXA2 ([RC211408])





Immunofluorescent staining of HepG2 cells using anti-FOXA2 mouse monoclonal antibody (TA500073) at 1:50 dilution.



HEK293T cells transfected with either pCMV6-ENTRY FOXA2 ([RC211408]) (Red) or empty vector control plasmid (Blue) were immunostained with anti-FOXA2 mouse monoclonal (TA500073, Dilution 1:1,000), and then analyzed by flow cytometry.