

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500121

CD31 (PECAM1) Mouse Monoclonal Antibody [Clone ID: OTI8A9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8A9
Applications: IHC, WB

Recommend Dilution: WB: 1:500, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 425-591 of human

PECAM1

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 82.4 kDa

Gene Name: platelet and endothelial cell adhesion molecule 1

Database Link: NP 000433 Entrez Gene 5175 Human

Background: The protein encoded by this gene is a cell adhesion molecule expressed on platelets and at

endothelial cellintercellular junctions

Synonyms: CD31; EndoCAM; endoCAM; GPIIA; PECA1; PECAM-1


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs), Leukocyte transendothelial migration






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PECAM1 ([RC208654], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PECAM1.



Immunohistochemical staining of paraffinembedded lung within the normal limits using anti-PECAM1 ([TA000121]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).