

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500245

FOXA1 Mouse Monoclonal Antibody [Clone ID: OTI3A8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3A8
Applications: IF, WB

Reactivity: WB 1:1000 Human, Dog

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 1-473 of human

FOXA1

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 49 kDa

Gene Name: forkhead box A1

Database Link: NP 004487 Entrez Gene 490657 DogEntrez Gene 3169 Human

Background: FOXA1 is a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear

factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver


Synonyms: HNF3A; TCF3A

Protein Families: ES Cell Differentiation/IPS, Transcription Factors







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FOXA1 ([RC206045], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FOXA1. Positive lysates [LY401427] (100ug) and [LC401427] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-FOXA1 monoclonal antibody.



Anti-FOXA1 mouse monoclonal antibody (TA500245) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FOXA1 ([RC206045]).