

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA500408**

## STK3 Mouse Monoclonal Antibody [Clone ID: OTI4G10]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI4G10
Applications: IHC, WB

**Reactivity:** WB 1:2000; IHC 1:50 Human, Monkey, Dog

Host: Mouse

Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human STK3 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 56.3 kDa

**Gene Name:** serine/threonine kinase 3

Database Link: NP 006272 Entrez Gene 475045 DogEntrez Gene 707425 MonkeyEntrez Gene 6788 Human

**Background:** Protein kinase activation is a frequent response of cells to treatment with growth

factors, chemicals, heat shock, or apoptosis-inducing agents. This protein kinase activation presumably allows cells to resist unfavorable environmental conditions. The yeast 'sterile 20'(Ste20) kinase acts upstream of the mitogen-activated protein kinase (MAPK) cascade that is activated under a variety of stress conditions. MST2 was identified as a kinase that is

activated by the proapoptotic agents straurosporine and FAS ligand

Synonyms: KRS1; MST2


**Protein Families:** Druggable Genome, Protein Kinase

**Protein Pathways:** MAPK signaling pathway







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY STK3 ([RC203897], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-STK3. Positive lysates [LY401892] (100ug) and [LC401892] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-STK3 monoclonal antibody.




Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-STK3 monoclonal antibody (1:200).






Western blot analysis of extracts (10ug) from 10 Human tissue by using anti-STK3 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon;10: spleen).




Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-STK3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of kidney tissue using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)




Immunohistochemical staining of paraffinembedded liver tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)






Immunohistochemical staining of paraffinembedded lung tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)



Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)





Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-STK3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500408, Dilution 1:50)