

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500515**

## LAT2 (SLC7A8) Mouse Monoclonal Antibody [Clone ID: OTI7D6]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI7D6
Applications: IHC, WB

**Recommend Dilution:** WB: 1:200 - 1:1000, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human SLC7A8 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.8 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 58.4 kDa

**Gene Name:** solute carrier family 7 member 8

Database Link: NP 036376 Entrez Gene 23428 Human

**Background:** Sodium-independent, high-affinity transport of small and large neutral amino acids such as

alanine, serine, threonine, cysteine, phenylalanine, tyrosine, leucine, arginine and tryptophan, when associated with SLC3A2/4F2hc. Acts as an amino acid exchanger. Has higher affinity for

L-phenylalanine than LAT1 but lower affinity for glutamine and serine. L-alanine is

transported at physiological concentrations. Plays a role in basolateral (re)absorption of neutral amino acids. Involved in the uptake of methylmercury (MeHg) when administered as

the L-cysteine or D,L-homocysteine complexes, and hence plays a role in metal ion homeostasis and toxicity. Involved in the cellular activity of small molecular weight nitrosothiols, via the stereoselective transport of L-nitrosocysteine (L-CNSO) across the transmembrane. Plays an essential role in the reabsorption of neutral amino acids from the

epithelial cells to the bloodstream in the kidney.


Synonyms: LAT2; LPI-PC1






**Protein Families:** Druggable Genome, Transmembrane


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SLC7A8 ([RC208586], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SLC7A8. Positive lysates [LY402177] (100ug) and [LC402177] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-SLC7A8 monoclonal antibody (1:200).




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-SLC7A8 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)






Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-SLC7A8mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-SLC7A8mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)




Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-SLC7A8mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)






Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of endometrium tissue using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)



Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-SLC7A8mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)





Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-SLC7A8mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500515, Dilution 1:50)