

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA500610**

Isocitrate dehydrogenase (IDH1) Mouse Monoclonal Antibody [Clone ID: OTI2H9]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2H9
Applications: IF, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human IDH1 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 46.7 kDa

**Gene Name:** isocitrate dehydrogenase (NADP(+)) 1, cytosolic

Database Link: NP 005887 Entrez Gene 3417 Human

**Background:** Isocitrate dehydrogenases catalyze the oxidative decarboxylation of isocitrate to 2-

oxoglutarate. These enzymes belong to two distinct subclasses, one of which utilizes NAD(+) as the electron acceptor and the other NADP(+). Five isocitrate dehydrogenases have been

reported: three NAD(+)-dependent isocitrate dehydrogenases, which localize to the

mitochondrial matrix, and two NADP(+)-dependent isocitrate dehydrogenases, one of which is mitochondrial and the other predominantly cytosolic. Each NADP(+)-dependent isozyme is

a homodimer. The protein encoded by this gene is the NADP(+)-dependent isocitrate

dehydrogenase found in the cytoplasm and peroxisomes. It contains the PTS-1 peroxisomal targeting signal sequence. The presence of this enzyme in peroxisomes suggests roles in the regeneration of NADPH for intraperoxisomal reductions, such as the conversion of 2, 4-

dienoyl-CoAs to 3-enoyl-CoAs, as well as in peroxisomal reactions that consume 2-oxoglutarate, namely the alpha-hydroxylation of phytanic acid. The cytoplasmic enzyme


serves a significant role in cytoplasmic NADPH production.




Synonyms: HEL-216; HEL-S-26; IDCD; IDH; IDP; IDPC; PICD

**Protein Pathways:** Citrate cycle (TCA cycle), Glutathione metabolism, Metabolic pathways


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IDH1 ([RC210582], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IDH1. Positive lysates [LY401782] (100ug) and [LC401782] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-IDH1 monoclonal antibody.




Western blot analysis of extracts (10ug) from 9 Human tissue by using anti-IDH1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: Colon).






Figure from citation: Western Blot of IDH1 protein level by using anti-IDH1 antibody in normal human astrocyte (NHA) cells. <u>View Citation</u>



Equivalent amounts of cell lysates (10 ug per lane) of wild-type Hela cells (WT, Cat# LC810HELA) and IDH1-Knockout Hela cells (KO, Cat# [LC810112]) were separated by SDS-PAGE and immunoblotted with anti-IDH1 monoclonal antibody TA500610. Then the blotted membrane was stripped and reprobed with anti-HSP90AB1 antibody ([TA500494]) as a loading control (1:500).



Anti-IDH1 mouse monoclonal antibody (TA500610) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IDH1 ([RC210582]).