

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA500618**

## SSX2 Mouse Monoclonal Antibody [Clone ID: OTI4A11]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4A11
Applications: IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human SSX2 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 21.6 kDa

**Gene Name:** SSX family member 2

Database Link: NP 003138 Entrez Gene 6757 Human

**Background:** The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX)

breakpoint proteins. These proteins may function as transcriptional repressors. They are also

capable of eliciting spontaneously humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. SSX1, SSX2 and SSX4 genes have been involved in the t(X;18) translocation characteristically found in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma translocation gene on chromosome 18 to one of the SSX genes on chromosome X. The encoded hybrid proteins are probably responsible for transforming activity. Two transcript

variants encoding distinct isoforms have been identified for this gene.


Synonyms: CT5.2; CT5.2A; HD21; HOM-MEL-40; SSX

Protein Families: Druggable Genome, Transcription Factors







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SSX2 ([RC201214], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SSX2. Positive lysates [LY418870] (100ug) and [LC418870] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-SSX2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)




Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)






Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)




Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500618, Dilution 1:50)





Anti-SSX2 mouse monoclonal antibody (TA500618) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY SSX2 ([RC201214]).