

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500757

ID3 Mouse Monoclonal Antibody [Clone ID: OTI15E3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI15E3

Applications: IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ID3 (NP_002158) produced in HEK293T

cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.64 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 13.0 kDa

Gene Name: inhibitor of DNA binding 3, HLH protein

Database Link: NP 002158 Entrez Gene 3399 Human

Background: Members of the ID family of helix-loop-helix (HLH) proteins lack a basic DNA-binding domain

and inhibit transcription through formation of nonfunctional dimers that are incapable of

binding to DNA

Synonyms: bHLHb25; HEIR-1


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors

Protein Pathways: TGF-beta signaling pathway







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ID3 ([RC200583], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ID3. Positive lysates [LY419492] (100ug) and [LC419492] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded breast tissue within the normal limits using anti-ID3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-ID3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)






Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)






Immunohistochemical staining of paraffinembedded liver tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)






Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)






Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-ID3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-ID3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500757, Dilution 1:50)





Anti-ID3 mouse monoclonal antibody (TA500757) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ID3 ([RC200583]).