

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501049

BHMT Mouse Monoclonal Antibody [Clone ID: OTI2F6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F6
Applications: IHC, WB

Recommend Dilution: WB 1:1000, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG3

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BHMT (NP_001704) produced in HEK293T

cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 44.8 kDa

Gene Name: betaine--homocysteine S-methyltransferase

Database Link: NP 001704 Entrez Gene 635 Human

Background: This gene encodes a cytosolic enzyme that catalyzes the conversion of betaine and

homocysteine to dimethylglycine and methionine, respectively. Defects in this gene could

lead to hyperhomocyst(e)inemia, but such a defect has not yet been observed.

Synonyms: BHMT1; HEL-S-61p


Protein Pathways: Cysteine and methionine metabolism, Glycine, serine and threonine metabolism, Metabolic

pathways







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BHMT ([RC203148], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BHMT. Positive lysates [LY400644] (100ug) and [LC400644] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-BHMT mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501049)



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-BHMT mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501049)





Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-BHMT mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501049)