

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501257

H6PD Mouse Monoclonal Antibody [Clone ID: OTI2A7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A7

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:250~500, IHC 1:50, IF 1:100, FLOW 1:100

Reactivity: Human, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human H6PD (NP_004276) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 88.7 kDa

Gene Name: hexose-6-phosphate dehydrogenase/glucose 1-dehydrogenase

Database Link: NP 004276 Entrez Gene 710107 MonkeyEntrez Gene 9563 Human

Background: There are 2 forms of glucose-6-phosphate dehydrogenase. G form is X-linked and H form,

encoded by this gene, is autosomally linked. This H form shows activity with other hexose-6-phosphates, especially galactose-6-phosphate, whereas the G form is specific for glucose-6-phosphate. Both forms are present in most tissues, but H form is not found in red cells.


Synonyms: CORTRD1; G6PDH; GDH

Protein Pathways: Metabolic pathways, Pentose phosphate pathway







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY H6PD ([RC209890], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-H6PD. Positive lysates [LY401369] (100ug) and [LC401369] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-H6PD monoclonal antibody.




Western blot analysis of extracts (10ug) from 10 Human tissue by using anti-H6PD monoclonal antibody at 1:500 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon;10: spleen).






Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-H6PD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501257, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-H6PD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501257, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-H6PD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501257, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-H6PD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501257, Dilution 1:50)




Anti-H6PD mouse monoclonal antibody (TA501257) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY H6PD ([RC209890]).



HEK293T cells transfected with either [RC209890] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-H6PD antibody (TA501257), and then analyzed by flow cytometry.





Flow cytometric Analysis of Jurkat cells, using anti-H6PD antibody (TA501257), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).