

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501353

CYB5R3 Mouse Monoclonal Antibody [Clone ID: OTI2A10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A10

Applications: FC, IHC, WB

Recommend Dilution: WB 1:200~500, IHC 1:150, FLOW 1:100

Reactivity: Human, Monkey, Dog

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CYB5R3 (NP_000389) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.76 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 34.1 kDa

Gene Name: cytochrome b5 reductase 3

Database Link: NP 000389 Entrez Gene 474479 DogEntrez Gene 714058 MonkeyEntrez Gene 1727 Human

Background: This gene encodes cytochrome b5 reductase, which includes a membrane-bound form in

somatic cells (anchored in the endoplasmic reticulum, mitochondrial and other membranes)

and a soluble form in erythrocytes. The membrane-bound form exists mainly on the

cytoplasmic side of the endoplasmic reticulum and functions in desaturation and elongation of fatty acids, in cholesterol biosynthesis, and in drug metabolism. The erythrocyte form is located in a soluble fraction of circulating erythrocytes and is involved in methemoglobin reduction. The membrane-bound form has both membrane-binding and catalytic domains, while the soluble form has only the catalytic domain. Alternate splicing results in multiple transcript variants. Mutations in this gene cause methemoglobinemias. [provided by RefSeq]

Synonyms: B5R; DIA1


Protein Families: Druggable Genome






Protein Pathways: Amino sugar and nucleotide sugar metabolism


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CYB5R3 ([RC201592], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CYB5R3. Positive lysates [LY400140] (100ug) and [LC400140] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CYB5R3 monoclonal antibody.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)




Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)






Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)






Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)





Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-CYB5R3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501353)









Flow cytometric Analysis of Jurkat cells, using anti-CYB5R3 antibody (TA501353), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).