

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA501363

MOB4A (MOB1B) Mouse Monoclonal Antibody [Clone ID: OTI3G9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3G9
Applications: IHC, WB

Recommend Dilution: WB 1:500, IHC 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MOBKL1A(NP_775739) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 24.9 kDa

Gene Name: MOB kinase activator 1B

Database Link: NP 775739 Entrez Gene 92597 Human

Background: The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds

Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint


regulation. [provided by RefSeq]. COMPLETENESS: complete on the 3' end.

Synonyms: MATS2; MOB4A; MOBKL1A






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MOBKL1A ([RC206337], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MOBKL1A. Positive lysates [LY406621] (100ug) and [LC406621] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-MOBKL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501363, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-MOBKL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501363, Dilution 1:50)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-MOBKL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501363, Dilution 1:50)



Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-MOBKL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501363, Dilution 1:50)





Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-MOBKL1A monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).