

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA502219**

# NMNAT1 Mouse Monoclonal Antibody [Clone ID: OTI1F7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1F7
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NMNAT1 (NP\_073624) produced in

HEK293T cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 31.8 kDa

**Gene Name:** nicotinamide nucleotide adenylyltransferase 1

Database Link: NP 073624 Entrez Gene 64802 Human

**Background:** The coenzyme NAD and its derivatives are involved in hundreds of metabolic redox reactions

and are utilized in protein ADP-ribosylation, histone deacetylation, and in some Ca(2+) signaling pathways. NMNAT (EC 2.7.7.1) is a central enzyme in NAD biosynthesis, catalyzing the condensation of nicotinamide mononucleotide (NMN) or nicotinic acid mononucleotide (NaMN) with the AMP moiety of ATP to form NAD or NaAD (Zhang et al., 2003 [PubMed

12574164]). [supplied by OMIM]


Synonyms: LCA9; NMNAT; PNAT1

**Protein Pathways:** Metabolic pathways, Nicotinate and nicotinamide metabolism






### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NMNAT1 ([RC204825], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NMNAT1. Positive lysates [LY402948] (100ug) and [LC402948] (20ug) can be purchased separately from OriGene.



Anti-NMNAT1 mouse monoclonal antibody (TA502219) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NMNAT1 ([RC204825]).