

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502263

DNAJB2 Mouse Monoclonal Antibody [Clone ID: OTI5E5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5E5
Applications: FC, WB

Recommend Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DNAJB2 (NP_001034639) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.48 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 30.4 kDa

Gene Name: Dnal heat shock protein family (Hsp40) member B2

Database Link: NP 001034639 Entrez Gene 3300 Human

Background: This gene is almost exclusively expressed in the brain, mainly in the neuronal layers. It

encodes a protein that shows sequence similarity to bacterial DnaJ protein and the yeast homologs. In bacteria, this protein is implicated in protein folding and protein complex dissociation. Alternatively spliced transcript variants encoding different isoforms have been


described for this gene. [provided by RefSeq, Jul 2011]

Synonyms: CMT2T; DSMA5; HSJ-1; HSPF3






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DNAJB2 ([RC206769], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNAJB2. Positive lysates [LY422070] (100ug) and [LC422070] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC206769] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DNAJB2 antibody (TA502263), and then analyzed by flow cytometry.