

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502413

Arginase 1 (ARG1) Mouse Monoclonal Antibody [Clone ID: OTI4E6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4E6

Applications: FC, IHC, WB

Recommend Dilution: WB 1:1000~2000, IHC 1:50, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ARG1 (NP_000036) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 34.6 kDa **Gene Name:** arginase 1

Database Link: NP 000036 Entrez Gene 383 Human

Background: Arginase catalyzes the hydrolysis of arginine to ornithine and urea. At least two isoforms of

mammalian arginase exist (types I and II) which differ in their tissue distribution, subcellular localization, immunologic crossreactivity and physiologic function. The type I isoform encoded by this gene, is a cytosolic enzyme and expressed predominantly in the liver as a component of the urea cycle. Inherited deficiency of this enzyme results in argininemia, an autosomal recessive disorder characterized by hyperammonemia. [provided by RefSeq]

Synonyms: arginase; arginase 1; liver; liver-type arginase; OTTHUMP00000017209; type I arginase


Protein Families: Druggable Genome

Protein Pathways: Arginine and proline metabolism, Metabolic pathways







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ARG1 ([RC204649], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ARG1. Positive lysates [LY424951] (100ug) and [LC424951] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ARG1 monoclonal antibody.




Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-ARG1 monoclonal antibody (1:200).





Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ARG1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502413)



HEK293T cells transfected with either [RC204649] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ARG1 antibody (TA502413), and then analyzed by flow cytometry.