

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA502521**

## NNMT Mouse Monoclonal Antibody [Clone ID: OTI3E3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3E3

**Applications:** FC, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NNMT (NP\_006160) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.49 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 29.4 kDa

Gene Name: nicotinamide N-methyltransferase

Database Link: NP 006160 Entrez Gene 4837 Human

**Background:** N-methylation is one method by which drug and other xenobiotic compounds are

metabolized by the liver. This gene encodes the protein responsible for this enzymatic activity

which uses S-adenosyl methionine as the methyl donor. [provided by RefSeq]


**Synonyms:** nicotinamide N-methyltransferase

**Protein Pathways:** Metabolic pathways, Nicotinate and nicotinamide metabolism







### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NNMT ([RC200641], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NNMT. Positive lysates [LY401860] (100ug) and [LC401860] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-NNMT mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502521)




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-NNMT mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502521)






HEK293T cells transfected with either [RC200641] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-NNMT antibody (TA502521), and then analyzed by flow cytometry.




Flow cytometric Analysis of Jurkat cells, using anti-NNMT antibody (TA502521), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).



Flow cytometric Analysis of Hela cells, using anti-NNMT antibody (TA502521), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).





Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-NNMT monoclonal antibody (1:200).